

# O-BRUSH

Surgical Masks ASTM Level 2

# TECHNICAL DATA SHEET



## Description

ASTM F2100-19 - Level 2

Bacterial Filtration Efficiency (BFE)  $\geq 98\%$

Particulate Filtration Efficiency at 0.1 micron (PFE)  $\geq 98\%$

Differential Pressure  $< 6\text{mm H}_2\text{O}/\text{cm}^2$

Resistance to penetration by synthetic blood – 120mm Hg

Flame Spread – Class 1

Latex Free

Single Use Only

Manufactured in China

## Dimensions

Length: 17.5 cm

Width: 9.5 cm

## Directions

- Clean your hands with soap and water or hand sanitizer before touching the mask.
- Remove a mask from the box and make sure there are no obvious tears or holes in either side of the mask.
- Determine which side of the mask is the top. The side of the mask that has a stiff bendable edge is the top and is meant to mold to the shape of your nose.
- The colored side of the mask is the front and should face away from you, while the white side touches your face.
- Hold the mask by the ear loops. Place a loop around each ear.
- Mold or pinch the stiff edge to the shape of your nose.
- Pull the bottom edge over the mask to cover your chin.

## Sterilisation

These products are non-sterile.

These products cannot be sterilized.

## Packaging

Shipping carton of 2000 masks.

50 masks per box and 40 boxes are placed within one carton box

Box Dimensions/Weight

52x39x30cm/ 7.5kg per carton

## Regulatory Information

Product CE marked as per 93/42/EEC Directive on Medical Devices.

Class of the device: I.

## Storage Information

Store in a dry and cool place, away from intense sources of heat. Keep the masks inside the box.

## Shelf Life

5 years

## Caution

When properly worn, this mask helps reduce exposure to blood and body fluids and minimize patient contamination. This product is not a respirator and does not eliminate exposure to or the risk of contracting disease or infection. Dispose if damaged or contaminated with blood or body fluids.

The dimensions and properties listed above can vary within pre-established specifications. This document was created using the most recent information. In the interest of continuous improvement, the characteristics of the product may change without prior notice.

ASTM F2100 STANDARD ASTM F2100 STANDARD ESTANDARD ASTM F2100			
ASTM F2100 Test Standard	Level 1	Level 2	Level 3
Particle Filtration Efficiency (PFE) <i>Efficacité de filtration des particules (EFP)</i> <i>Eficiencia de filtración de partículas</i>	≥95%	≥98%	≥98%
Bacterial Filtration Efficiency (BFE) <i>Efficacité de filtration des bactéries (EFB)</i> <i>Eficiencia de filtración de bacterias</i>	≥95%	≥98%	≥98%
ΔP (mmHg/cm²)	<4.0	<5.0	<5.0
Flame Spread <i>Classement d'inflammabilité</i> <i>Clasificación de inflamabilidad</i>	Class 1	Class 1	Class 1
Resistance to penetration by synthetic blood <i>Résistance à la pénétration du sang synthétique</i> <i>Resistencia a la penetración de sangre sintética</i>	80 mmHg	120 mmHg	160 mmHg



## Bacterial Filtration Efficiency (BFE) and Differential Pressure (Delta P) Final Report

Test Article: Surgical Face Mask, FM-34EE/MK7200  
LOT #191028  
Purchase Order: [REDACTED]  
Study Number: [REDACTED]  
Study Received Date: 04 Dec 2019  
Testing Facility: Nelson Laboratories, LLC  
6280 S. Redwood Rd.  
Salt Lake City, UT 84123 U.S.A.  
Test Procedure(s): Standard Test Protocol (STP) Number: STP0004 Rev 18  
Deviation(s): None


**Summary:** The BFE test is performed to determine the filtration efficiency of test articles by comparing the bacterial control counts upstream of the test article to the bacterial counts downstream. A suspension of *Staphylococcus aureus* was aerosolized using a nebulizer and delivered to the test article at a constant flow rate and fixed air pressure. The challenge delivery was maintained at  $1.7 - 3.0 \times 10^3$  colony forming units (CFU) with a mean particle size (MPS) of  $3.0 \pm 0.3 \mu\text{m}$ . The aerosols were drawn through a six-stage, viable particle, Andersen sampler for collection. This test method complies with ASTM F2101-19 and EN 14683:2019, Annex B.


The Delta P test is performed to determine the breathability of test articles by measuring the differential air pressure on either side of the test article using a manometer, at a constant flow rate. The Delta P test complies with EN 14683:2019, Annex C and ASTM F2100-19.

All test method acceptance criteria were met. Testing was performed in compliance with US FDA good manufacturing practice (GMP) regulations 21 CFR Parts 210, 211 and 820.

Test Side: Inside  
BFE Test Area:  $\sim 40 \text{ cm}^2$   
BFE Flow Rate: 28.3 Liters per minute (L/min)  
Delta P Flow Rate: 8 L/min  
Conditioning Parameters:  $85 \pm 5\%$  relative humidity (RH) and  $21 \pm 5^\circ\text{C}$  for a minimum of 4 hours  
Test Article Dimensions:  $\sim 173 \text{ mm} \times \sim 172 \text{ mm}$   
Positive Control Average:  $1.9 \times 10^3 \text{ CFU}$   
Negative Monitor Count:  $< 1 \text{ CFU}$   
MPS:  $3.0 \mu\text{m}$



  
Study Director

  
Janelle R. Bentz, M.S.

18 Dec 2019  
Study Completion Date

**Results:**

Test Article Number	Percent BFE (%)
1	99.7
2	99.6
3	99.8
4	99.8
5	>99.9

Test Article Number	Delta P (mm H <sub>2</sub> O/cm <sup>2</sup> )
1	4.6
2	4.4
3	5.2
4	5.0
5	4.7

The filtration efficiency percentages were calculated using the following equation:

$$\% BFE = \frac{C - T}{C} \times 100$$

C = Positive control average

T = Plate count total recovered downstream of the test article

Note: The plate count total is available upon request

## Latex Particle Challenge GLP Report

Test Article: Surgical Face Mask, FM-34EE  
LOT #191028  
Purchase Order: [REDACTED]  
Study Number: [REDACTED]  
Study Received Date: 05 Nov 2019  
Testing Facility: Nelson Laboratories, LLC  
6280 S. Redwood Rd.  
Salt Lake City, UT 84123 U.S.A.  
Test Procedure(s): Standard Test Protocol (STP) Number: STP0005 Rev 07  
Deviation(s): None

**Summary:** This procedure was performed to evaluate the non-viable particle filtration efficiency (PFE) of the test article. Monodispersed polystyrene latex spheres (PSL) were nebulized (atomized), dried, and passed through the test article. The particles that passed through the test article were enumerated using a laser particle counter.

Three one-minute counts were performed, with the test article in the system, and the results averaged. Three one-minute control counts were performed, without a test article in the system, before and after each test article and the counts were averaged. Control counts were performed to determine the average number of particles delivered to the test article. The filtration efficiency was calculated using the average number of particles penetrating the test article compared to the average of the control values.

The procedure employed the basic particle filtration method described in ASTM F2299, with some exceptions; notably the procedure incorporated a non-neutralized challenge. In real use, particles carry a charge, thus this challenge represents a more natural state. The non-neutralized aerosol is also specified in the FDA guidance document on surgical face masks. All test method acceptance criteria were met.

Test Side: Inside  
Area Tested: 91.5 cm<sup>2</sup>  
Particle Size: 0.1 µm  
Laboratory Conditions: 21°C, 23% relative humidity (RH) at 0833; 21°C, 24% RH at 1045  
Average Filtration Efficiency: 99.60%  
Standard Deviation: 0.270




*Curtis Gerow*  
Study Director

Curtis Gerow, B.S.

*13 Nov 2019*  
Study Completion Date





**Results:**

Test Article Number	Average Test Article Counts	Average Control Counts	Filtration Efficiency (%)
1	37	11,434	99.68
2	33	11,749	99.72
3	27	12,991	99.79
4	42	13,815	99.70
5	122	13,994	99.13

**Test Method Acceptance Criteria:** Ambient background particles detected through the test system must be below 1% of the challenge total (<100 particles).

**Procedures:**

**Test Set-up:** Testing was conducted in an ISO Class 5 (class 100) HEPA filtered hood. The inlet air to the test system was filtered through a 0.2 µm rated air filter. The particle generator outlet was clamped off and the number of background particles within the test system was verified to be <100 particles at 1 cubic foot per minute (CFM). The flow rate through the test system was maintained at 1 CFM ± 5%.

An aliquot of the PSL was aerosolized using a particle generator, mixed with additional filtered air, dried and passed through the test system. The particles delivered were enumerated using a laser based particle counter.

**Test Procedure:** A test article was placed into the holder and the system was allowed to stabilize. The average number of particles being delivered to the test article was determined (no medium in air stream) as triplicate one-minute control readings were taken prior to and after every test article. Control count averages were maintained at a level of 10,000-15,000 particles per cubic foot. Triplicate one-minute counts were recorded for the test article between the control counts.

The PFE of each test article was determined by using the following equation:

$$\% PFE = \frac{C - T}{C} \times 100$$

Where: C = Combined average of the control counts  
T = Average test article counts

## Quality Assurance Statement

**Compliance Statement:** The test was conducted in accordance with the USFDA (21 CFR Parts 58, 210, 211, and 820) Regulations. This final report reflects the raw data.

Activity	Date
Study Initiation	11 Nov 2019
Phase Inspected by Quality Assurance: Latex Test	12 Nov 2019
Audit Results Reported to Study Director	12 Nov 2019
Audit Results Reported to Management	12 Nov 2019

Scientists	Title
Sarah Smit	Supervisor
Curtis Gerow	Study Director

**Data Disposition:** The study plan, raw data and final report from this study are archived at Nelson Laboratories, LLC or an approved off-site location.

  
\_\_\_\_\_  
Quality Assurance

15 Nov 2019  
\_\_\_\_\_  
Date

## Synthetic Blood Penetration Resistance Final Report

Test Article: Surgical Face Mask, FM-34EE/MK7200  
LOT #191028  
Purchase Order: [REDACTED]  
Study Number: [REDACTED]  
Study Received Date: [REDACTED]  
Testing Facility: Nelson Laboratories, LLC  
6280 S. Redwood Rd.  
Salt Lake City, UT 84123 U.S.A.  
Test Procedure(s): Standard Test Protocol (STP) Number: STP0012 Rev 09  
Deviation(s): None

**Summary:** This procedure was performed to evaluate surgical facemasks and other types of protective clothing materials designed to protect against fluid penetration. The purpose of this procedure is to simulate an arterial spray and evaluate the effectiveness of the test article in protecting the user from possible exposure to blood and other body fluids. The distance from the target area surface to the tip of the cannula is 30.5 cm. A test volume of 2 mL of synthetic blood was employed using the targeting plate method.

This test method was designed to comply with ASTM F1862 and ISO 22609 (as referenced in EN 14683:2019 and AS4381:2015) with the following exception: ISO 22609 requires testing to be performed in an environment with a temperature of  $21 \pm 5^{\circ}\text{C}$  and a relative humidity of  $85 \pm 10\%$ . Instead, testing was performed at ambient conditions within one minute of removal from the environmental chamber held at those parameters.

All test method acceptance criteria were met. Testing was performed in compliance with US FDA good manufacturing practice (GMP) regulations 21 CFR Parts 210, 211 and 820.


Number of Test Articles Tested: 32  
Number of Test Articles Passed: 30  
Test Side: Outside  
Pre-Conditioning: Minimum of 4 hours at  $21 \pm 5^{\circ}\text{C}$  and  $85 \pm 5\%$  relative humidity (RH)  
Test Conditions:  $21.8^{\circ}\text{C}$  and 22% RH

**Results:** Per ASTM F1862 and ISO 22609, an acceptable quality limit of 4.0% is met for a normal single sampling plan when  $\geq 29$  of 32 test articles show passing results.

Test Pressure: 120 mmHg (16.0 kPa)

Test Article Number	Synthetic Blood Penetration
1-7, 9-27, 29-32	None Seen
8, 28	Yes

  
Study Director

  
For  
Janelle R. Bentz, M.S.

  
Study Completion Date